# M5 Super LipoRNAiMAX Transfection Reagent 使用说明书

产品名称	单位	货号
M5 Super LipoRNAiMAX Transfection Reagent	0.25 ml	MF139-01
M5 Super LipoRNAiMAX Transfection Reagent	0.50 ml	MF139-02

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**B、 Transfection of siRNA**

Reverse transfection is faster to perform than forward transfection and is the method of choice for high-throughput transfection. Perform reverse transfection by preparing complexes inside the wells, and then adding cells and medium. Because the cells and siRNA-lipid complexes are prepared on the same day, we recommended using 2.5 $\times$  more cells than for a regular transfection method.

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Transfect cells according to the following table. The transfection is designed for one RNA amount combined with one amount of LipoRNAiMAX. The prepared mix is enough to have triplicates (96-well), duplicates (24-well), and single well (6-well) transfections, and account for pipetting variations. For additional information on scaling your transfection reaction.

Timeline		Steps	Procedure Details			
Day 0	1	Seed cells to be 60-80% confluent at transfection	Component	96-well	24-well	6-well
	2	Dilute LipoRNAiMAX Reagent in Opti-MEM <sup>®</sup> Medium	Adherent cells	1-4 $\times 10^4$	0.5-2 $\times 10^5$	0.25-1 $\times 10^6$
	3	Dilute siRNA in Opti-MEM <sup>®</sup> Medium	Opti-MEM <sup>®</sup> Medium	25 $\mu$ L	50 $\mu$ L	150 $\mu$ L
Day 1	4	Add diluted siRNA to diluted LipoRNAiMAX Reagent (1:1 ratio)	Lipo RNAiMAX Reagent	1.5 $\mu$ L	3 $\mu$ L	9 $\mu$ L
	5	Incubate	Opti-MEM <sup>®</sup> Medium	25 $\mu$ L	50 $\mu$ L	150 $\mu$ L
	6	Add siRNA-lipid complex to cells	siRNA (10 $\mu$ M)	0.5 $\mu$ L (5 pmol)	1 $\mu$ L (10 pmol)	3 $\mu$ L (30 pmol)
	7	Visualize/analyze transfected cells	Diluted siRNA	25 $\mu$ L	50 $\mu$ L	150 $\mu$ L
Day 2-4			Diluted LipoRNAiMAX Reagent	25 $\mu$ L	50 $\mu$ L	150 $\mu$ L
			Incubate for 5 minutes at room temperature.			
			Component	96-well	24-well	6-well
			siRNA-lipid complex per well	10 $\mu$ L	50 $\mu$ L	250 $\mu$ L
		Final siRNA used per well	1 pmol	5 pmol	25 pmol	
		Final LipoRNAiMAX used per well	0.3 $\mu$ L	1.5 $\mu$ L	7.5 $\mu$ L	
			Incubate cells for 1-3 days at 37°C. Then, analyze transfected cells.			

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