M5 Hiper Multi-color Self-Lighting Western Protein Marker (10-200 kDa)

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
M5 Hiper Multi-color Self-Lighting Western Protein Marker (10-200 kDa) is a ready-to-use mixture with 10 IgG-binding proteins covering a wide range of molecular weights from 10 to 200 kDa in Tris Glycine buffer. M5 Hiper Multi-color Self-Lighting Western Protein Marker (10-200 kDa) performs dual functions. First, it contains 11 pre-stained proteins (10, 15, 25, 35, 45, 65, **70, 75**, 100, 135 and 180 kDa) for monitoring protein separation during SDS-PAGE, verification of Western transfer efficiency on membranes (nitrocellulose, PVDF, or nylon) and for approximating the protein size. Second, 10 IgG binding proteins can be immuno-detected on film or by CCD imaging. M5 Hiper Multi-color Self-Lighting Western Protein Marker (10-200 kDa) is compatible for chemiluminescent, fluorescent, chromogenic or other detection systems. In addition, M5 Hiper Multi-color Self-Lighting Western Protein Marker (10-200 kDa) has two reference bands with enhanced intensity (at 30 kDa and 80 kDa). The marker is supplied in the gel loading buffer and is ready to use. Do NOT heat, dilute, or add reducing agents before loading.

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- a. Transfer with buffer containing 20% methanol to fix proteins on membrane;
- b. Wash membrane with PBS or TBS containing less than 0.1% Tween-20 at 4 $^{\circ}$ C.

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