



M5 Anti-DYKDDDDK-Tag mAb (Magnetic Beads)

Catalog Number MF098-plus

Storage Temperature $-20\text{ }^{\circ}\text{C}$

TECHNICAL BULLETIN

Product Description

ANTI-FLAG M2 Magnetic Beads are composed of the murine derived, ANTI-FLAG M2 monoclonal antibody attached to superparamagnetic iron impregnated, 4% agarose beads with an average diameter of 50 μm . The M2 antibody binds to fusion proteins containing the FLAG[®] peptide sequence.¹ Additionally, the M2 antibody recognizes the FLAG octapeptide sequence (N-Asp-Tyr-Lys-Asp-Asp-Asp-Lys-C) at the N-terminus, Met-N-terminus, or C-terminus locations of a fusion protein in mammalian and bacterial extracts.

The ANTI-FLAG M2 Magnetic Beads are useful for detection and capture of fusion proteins containing a FLAG peptide sequence by commonly used immunoprecipitation procedures. The magnetic properties allow for very rapid separation of the beads from a suspension, significantly accelerating manipulations, such as repetitive washings or processing of multiple samples performed in multiwell plates. This leads to faster experimentation, better reproducibility, and more accurate quantitation of the proteins of interest.

Binding capacity: 0.6 mg of FLAG fusion protein per 1 ml of packed magnetic beads.

Specificity: $\geq 90\%$ specificity towards FLAG fusion proteins from mammalian and bacterial cell extracts.

Reagent

The ANTI-FLAG M2 Magnetic Bead resin is supplied as a 50% suspension in 50% glycerol with 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.4, and 0.02% (w/v) sodium azide (PBA/A).

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Equipment Required But Not Provided

Magnetic Separators for:

Microcentrifuge tubes (Catalog Number M1167)

Tissue culture flasks (Catalog Number M1292)

Centrifuge tubes (Catalog Number M1542)

Magnet for 96-well tissue culture plates

(Catalog Number SHM05)

Magnetic plate for standard sized well plates, T-25

through T-75 tissue culture flasks, and up to 5 cm dishes (Catalog Number SHM04)

Do not use a magnetic stirring system. This will destroy the resin beads.

Storage/Stability

The ANTI-FLAG M2 Magnetic Beads ship on wet ice and storage at $-20\text{ }^{\circ}\text{C}$ is recommended. The product is supplied in a 50% glycerol solution with preservative and is stable for 2 years at $-20\text{ }^{\circ}\text{C}$. After use, the resin should be cleaned and stored in 50% glycerol with TBS or PBS buffer containing preservative to protect the product.

Freezing the magnetic beads in the absence of 50% glycerol will irreversibly damage the bead structure.

Procedures

Note: It is recommended the entire technical bulletin be read before use, especially the Reagent Compatibility Table at the end of this bulletin.

There are many different procedures for performing small-scale affinity capture experiments. The following procedures are written for a single sample.

For batch-wise purification, 100 μl of the resin bead suspension per reaction ($\sim 50\text{ }\mu\text{l}$ of packed gel) is recommended. For use in a 96-well plate format, 10 μl of packed gel is recommended per well. The amount of resin bead can be varied depending on the amount of target protein in the sample and the type of magnetic separator utilized.

Part I. Sample Preparation

1. Adjust the pH of the protein extract to between pH 7–8. It is also useful to adjust the salt concentration with sodium or potassium chloride to ≥ 0.15 M in the protein extract to prevent a large amount of nonspecific protein binding to the resin.
2. The FLAG fusion protein extract must be clarified to remove any insoluble material. A large amount of insoluble material may require centrifugation (10,000–20,000 $\times g$ for 15 minutes) for removal. The protein extract should also be filtered through a 0.45 or 0.22 μm filter to remove any remaining cell debris and particulates that may interfere with protein binding.

Part II. Binding Procedures

For purification of FLAG fusion proteins, the resin can be used in either a batch or 96-well plate format. For larger volumes of lysate, the batch format is recommended to quickly capture the target protein from a large volume of extract. If a smaller sample is being purified, the FLAG fusion protein can be immuno-precipitated.

A. Automation Format for ANTI-FLAG M2 Magnetic Beads - Please see our website (www.sigma.com/automation) for a purification procedure for FLAG fusion proteins using the KingFisher[®] automated platform.

B. Batch Format for Absorption of FLAG Fusion Proteins with ANTI-FLAG M2 Magnetic Beads - This procedure provides a quick and efficient way to purify FLAG fusion proteins from a dilute solution.

Table 1.
Binding Capacity of Beads

Packed Gel Volume (μl)	Binding Capacity (μg)
20	~12
50	~30
100	~60
200	~120

The ANTI-FLAG M2 Magnetic Beads are stored in 50% glycerol with buffer. The glycerol must be removed just prior to use and the resin equilibrated with buffer (steps 1–5). The equilibration can be done at room temperature or at 2–8 °C. Remove only the volume of resin that is necessary for the purification (see Table 1).

Do not allow the resin to remain in TBS (50 mM Tris HCl, 150 mM NaCl, pH 7.4) buffer for extended periods of time (>24 hours) unless an antimicrobial agent (e.g., 0.02% sodium azide) is added to the buffer.

1. Thoroughly resuspend the resin by gentle inversion. Make sure the bottle of ANTI-FLAG M2 Magnetic Beads is a uniform suspension. Remove an appropriate volume for use (see Table 1).
2. Transfer resin to an appropriate sized tube. Equilibrate beads by resuspending with 5 packed gel volumes of TBS (50 mM Tris HCl, 150 mM NaCl, pH 7.4). Mix thoroughly. Place tube in the appropriate magnetic separator to collect the beads. Remove and discard the storage buffer/TBS mixture.
3. Equilibrate beads by resuspending with another 5 packed gel volumes of TBS (50 mM Tris HCl, 150 mM NaCl, pH 7.4). Mix thoroughly.
4. Place tube in the appropriate magnetic separator to collect the beads. Remove and discard TBS buffer.
5. Repeat steps 3 and 4 once. Allow a small amount of buffer to remain on the top of the beads.
6. Incubate the protein extract (see Sample Preparation) with the equilibrated beads (step 5) for ~1 hour at room temperature with gentle mixing to capture the FLAG fusion proteins. Mixing should be done on either a rotating device or a platform shaker. **Do not use a magnetic stirring system.** This will destroy the resin beads.
7. Once the binding step is complete, collect the magnetic beads by placing the tube in the appropriate magnetic separator and remove the supernatant.
8. Wash the resin beads with TBS buffer to remove all of the nonspecifically bound proteins. Washing should be done with 20 packed gel volumes of TBS buffer, performed in three sequential bead washing. **Note:** The washing process can be monitored by measuring the absorbance of the supernatant at 280 nm. Continue washing the resin until the absorbance difference between the wash solution aspirated from the beads and the wash solution (TBS) blank is < 0.05 .

9. The FLAG proteins can be eluted from the magnetic beads either by low pH or by competition with the FLAG peptide.
- Elution of FLAG fusion proteins under acidic conditions with glycine - Elute the bound FLAG fusion protein from the magnetic beads with 10 packed gel volumes of 0.1 M glycine HCl, pH 3.0, collecting 1 packed gel volume of eluate in a vial containing 15–25 μ l of 1 M Tris, pH 8.0. Do not leave eluate in the glycine HCl solution for longer than 20 minutes. Re-equilibrate the resin to neutral pH as soon as possible after elution.
Or
 - Elution of FLAG fusion proteins by competition with the FLAG Peptide – Elute the bound FLAG fusion protein by competitive elution with five packed gel volumes of a solution containing the FLAG peptide (100 μ g/ml, Catalog Number F3290) in TBS buffer.
10. Cleaning the Magnetic Beads - It is recommended the beads be cleaned immediately after use by washing with three packed gel volumes of 0.1 M glycine HCl, pH 3.0. The beads should be immediately re-equilibrated in TBS buffer until the effluent is at neutral pH.
11. Storing the Magnetic Beads - After cleaning the beads, collect the magnetic beads by placing the tube in the appropriate magnetic separator and remove the buffer. The beads may be stored as a 50% suspension in 50% glycerol with TBS or PBS buffer containing 0.02% sodium azide. Store the beads at 2–8 °C or –20 °C without draining.

C. Format for Immunoprecipitation of FLAG Fusion Proteins using ANTI-FLAG M2 Magnetic Beads - This procedure is recommended for the purification of small amounts of FLAG fusion proteins.

Note: For antigens and protein:protein complexes requiring a special lysis buffer composed of a different percentage of detergent, it is recommended to pretest the resin before use. The ANTI-FLAG M2 resin bead is resistant to the many detergents at the following concentrations: 5.0% TWEEN® 20, 5.0% Triton™ X-100, 0.1% IGEPAL® CA-630, 0.1% CHAPS, and 0.2% digitonin. It can also be used with 1.0 M NaCl or 1.0 M urea. See the Reagent Compatibility Table for use with additional chemicals.

The following procedure is an example of a single immunoprecipitation reaction. For multiple immunoprecipitation reactions, calculate the volume of reagents needed according to the number of samples to be processed. For immunoprecipitation reactions, it is recommended to use 40 μ l of the 50% bead

suspension per reaction (~20 μ l of packed gel volume). Smaller amounts of resin (~10 μ l of packed gel volume, which binds >1 μ g FLAG fusion protein) can be used.

Note: Two control reactions are recommended for the procedure. The first control is immunoprecipitation with FLAG-BAP™ fusion protein (positive control) and the second is a reagent blank with no protein (negative control).

- Thoroughly resuspend the resin by gentle inversion. Make sure the bottle of ANTI-FLAG M2 Magnetic Beads is a uniform suspension. Remove an appropriate volume for use (see Table 1). To lessen damage to beads, it is recommended to cut the end of the pipette tip.
- Place tube in the appropriate magnetic separator to collect the beads. Aspirate and discard storage buffer.
- Wash the packed gel twice with 10 packed gel volumes of TBS (50 mM Tris HCl, 150 mM NaCl, pH 7.4) buffer. Be sure most of the wash buffer is removed and no resin is discarded.
Note: For multiple immunoprecipitation samples, wash the total volume of resin needed for all samples together. After washing, resuspend the resin in TBS buffer and divide the resin according to the number of samples tested. Place tube in the appropriate magnetic separator to collect the beads. Remove and discard TBS buffer.
- Add 200–1,000 μ l of cell lysate to the washed resin beads. If necessary, bring the final volume to 1 ml by adding lysis buffer (50 mM Tris HCl, pH 7.4, with 150 mM NaCl, 1 mM EDTA, and 1% TRITON X-100). The volume of cell lysate to be used depends on the expression level of FLAG fusion protein in the transfected cells. For the positive control, add 1 ml of TBS buffer and 4 μ l of 50 ng/ μ l FLAG-BAP fusion protein (~200 ng) to the washed resin beads. For the negative control, add 1 ml of lysis buffer only with **no** protein. The amount of FLAG-BAP fusion protein to be precipitated depends on the detection method. 200 ng of protein is sufficient for an activity assay or for an immunoblot analysis. For SDS-PAGE analysis with Coomassie® blue or silver staining detection, use 1 μ g of FLAG-BAP fusion protein.
- Agitate or shake (a roller shaker is recommended) all samples and controls gently for 2 hours. In order to increase the binding efficiency, the binding step may be extended overnight.
- Place tubes in the appropriate magnetic separator to collect the beads and remove the supernatant with a narrow-end pipette tip.

7. Wash the resin three times with a total of 20 packed gel volumes of TBS (50 mM Tris HCl, 150 mM NaCl, pH 7.4). Make sure all the supernatant is removed by using a Hamilton® syringe or equivalent device.
8. Elution of the FLAG fusion proteins - Three elution methods are recommended according to protein characteristics or further usage:
 - Protein elution under native conditions by competition with 3X FLAG peptide. The elution efficiency is very high using this method.
 - Elution under acidic conditions with 0.1 M glycine HCl, pH 3.0. This is a fast and efficient elution method. Neutralization of the eluted protein with wash buffer may help preserve its activity.
 - Elution with sample buffer for gel electrophoresis and immunoblotting.

Elution with 3X FLAG peptide

- a. Prepare 3X FLAG elution solution. Dissolve 3X FLAG peptide (Catalog Number F4799) in 0.5 M Tris HCl, pH 7.5, with 1 M NaCl at a concentration of 25 µg/µl. Dilute 5-fold with water to prepare a 3X FLAG stock solution containing 5 µg/µl of 3X FLAG peptide. For elution, add 3 µl of 5 µg/µl 3X FLAG peptide stock solution to 100 µl of TBS buffer (150 ng/µl final concentration).
- b. Add 5 packed gel volumes of 3X FLAG elution solution to each sample and control resin.
- c. Incubate the samples and controls with gently shaking or on a rotator for 30 minutes at 2–8 °C.
- d. Place tube in the appropriate magnetic separator to collect the beads. Transfer the supernatants to fresh tubes using a Hamilton syringe or equivalent device. Be careful not to transfer any resin.
- e. Repeat steps a-d, pooling eluates in the same tube.
- f. For immediate use, store the combined eluates at 2–8 °C. Store at –20 °C for long term storage.
- g. For cleaning and storage of used resin, see Batch Format Procedure, steps 10 and 11.

Elution with 0.1 M Glycine HCl, pH 3.0 - The procedure should be performed at room temperature. **Do not leave the resin in this buffer more than 20 minutes.**

- a. Add 5 packed gel volumes of 0.1 M glycine HCl buffer, pH 3.0, to each sample and control resin.

- b. Incubate the samples and controls with gentle shaking or on a rotator for 5 minutes at room temperature.
- c. Place tube in the appropriate magnetic separator to collect the beads. Transfer the supernatants to fresh tubes containing 10 µl of 0.5 M Tris HCl, pH 7.4, with 1.5 M NaCl, using a Hamilton syringe or equivalent device. Be careful not to transfer any resin.
- d. Repeat steps a-d, pooling eluates in same tube.
- e. For immediate use, store the combined eluates at 2–8 °C. Store at –20 °C for long term storage.
- f. For cleaning and storage of used resin, see Batch Format Procedure, steps 10 and 11.

Elution with SDS-PAGE Sample Buffer - The procedure should be performed at room temperature. Sample buffer should be at room temperature before use. In order to minimize the denaturation and elution of the M2 antibody, no reducing agents, e.g., 2-mercaptoethanol or DTT, should be included in the sample buffer. The addition of reducing agents will result in the dissociation of the heavy and light chains of the immobilized M2 antibody (25 and 50 kDa bands). If reducing conditions are absolutely necessary, a reducing agent may be added. The final concentration of 2-mercaptoethanol or DTT in the 1× sample buffer (62.5 mM Tris HCl, pH 6.8, 2% SDS, 10% (v/v) glycerol, and 0.002% bromphenol blue) should be 5% or 50 mM, respectively. **Note:** Elution of the bound FLAG fusion protein as a SDS-PAGE sample results in damage to the ANTI-FLAG M2 Magnetic Beads and they cannot be used again. The SDS in the sample buffer will denature the M2 antibody and boiling will damage the bead structure.

- a. Add 20 µl of 2× sample buffer (125 mM Tris HCl, pH 6.8, 4% SDS, 20% (v/v) glycerol, and 0.004% bromphenol blue) to each sample and control.
- b. Boil the sample and control tubes for 3 minutes.
- c. Place tubes in the appropriate magnetic separator to collect the beads. Transfer the supernatants to fresh tubes with a Hamilton syringe or a narrow-end Pasteur pipette. The samples and controls are ready for loading on SDS-PAGE and immunoblotting using ANTI-FLAG or specific antibodies against the fusion protein.

References

1. Brizzard, B.L., et al., Immunoaffinity purification of FLAG epitope-tagged bacterial alkaline phosphatase using a novel monoclonal antibody and peptide elution. *BioTechniques*, **16**, 730-735 (1994).
2. Safarik, I., and Safarikova, M., Magnetic Techniques for the Isolation and Purification of Proteins and Peptides. *BioMagnetic Research and Technology*, **2**, (2004).

Reagent Compatibility Table

Reagent	Effects	Comments
Chaotropic agents (e.g., guanidine HCl or urea)	Denatures the immobilized M2 antibody	Do not use any reagent that contains these types of components since it will denature the M2 antibody on the resin and destroy its ability to bind the FLAG fusion proteins. Low concentrations of urea (1 M or less) can be used.
Reducing agents (such as 2-mercaptoethanol, DTT, or DTE)	Reduces the disulfide bridges holding the M2 antibody chains together	Do not use any reagent that contains these types of components since it will reduce the disulfide linkages in the M2 antibody on the resin and destroy its ability to bind the FLAG fusion proteins.
TWEEN 20, 5% or less	Reduces non-specific protein binding to the resin bead	May be used up to recommended concentration of 5%, but do not exceed.
TRITON X-100, 5% or less	Reduces non-specific protein binding to the resin bead	May be used up to recommended concentration of 5%, but do not exceed.
IGEPAL CA-630, 0.1% or less	Reduces non-specific protein binding to the resin bead	May be used up to recommended concentration of 0.1%, but do not exceed.
CHAPS, 0.1% or less	Reduces non-specific protein binding to the resin bead	May be used up to recommended concentration of 0.1%, but do not exceed.
Digitonin, 0.2% or less	Reduces non-specific protein binding to the resin bead	May be used up to recommended concentration of 0.2%, but do not exceed.
Sodium chloride, 1.0 M or less	Reduces non-specific protein binding to the resin bead by reducing ionic interactions	May be used up to recommended concentration of 1.0 M, but do not exceed.
Sodium dodecyl sulfate (SDS)	Denatures the immobilized M2 antibody	Do not use any reagent that contains this detergent in the loading and washing buffers since it will denature the M2 antibody on the resin bead and destroy its ability to bind the FLAG fusion proteins. SDS in the sample buffer is useful for the removal of proteins for immunoprecipitation, but the resin bead cannot be reused
0.1 M glycine HCl, pH 3.5	Elutes FLAG protein from the resin bead	Do not leave the column in glycine HCl for longer than 20 minutes. Longer incubation times will begin to denature the M2 antibody
Deoxycholate	Interferes with M2 binding to FLAG fusion proteins	Do not use any reagent that contains this detergent since it will inhibit the M2 antibody from binding to FLAG fusion proteins.

Troubleshooting Guide

Problem	Possible Cause	Solution
No signal is observed.	FLAG fusion protein is not present in the sample.	<ul style="list-style-type: none"> • Make sure the protein of interest contains the FLAG sequence by immunoblot or dot blot analyses. • Prepare fresh lysates. Avoid using frozen lysates. • Use appropriate protease inhibitors in the lysate or increase their concentrations to prevent degradation of FLAG fusion protein
	Washes are too stringent.	<ul style="list-style-type: none"> • Reduce the number of washes. • Avoid adding high concentrations of NaCl to the mixture. • Use solutions that contain less or no detergent.
	Incubation times are inadequate.	<ul style="list-style-type: none"> • Increase the incubation times with the affinity resin (from several hours to overnight).
	Interfering substance is present in sample.	<ul style="list-style-type: none"> • Lysates containing high concentrations of dithiothreitol (DTT), 2-mercaptoethanol, or other reducing agents may destroy antibody function, and must be avoided. • Excessive detergent concentrations may interfere with the antibody-antigen interaction. Detergent levels in buffers may be reduced by dilution.
	Detection system is inadequate.	<p>If Western blotting detection is used:</p> <ul style="list-style-type: none"> • Check primary and secondary antibodies using proper controls to confirm binding and reactivity. • Verify the transfer was adequate by staining the membrane with Ponceau S. • Use fresh detection substrate or try a different detection system
Background is too high.	Proteins bind nonspecifically to the ANTI-FLAG monoclonal antibody, the resin beads, or the microcentrifuge tubes	<ul style="list-style-type: none"> • Pre-clear lysate with Mouse IgG-Agarose (Catalog Number A0919) to remove nonspecific binding proteins. • After suspending beads for the final wash, transfer entire sample to a clean microcentrifuge tube before centrifugation
	Washes are insufficient.	<ul style="list-style-type: none"> • Increase the number of washes. • Increase duration of the washes, incubating each wash for at least 15 minutes. • Increase the salt and/or detergent concentrations in the wash solutions. • Centrifuge at lower speed to avoid nonspecific trapping of denatured proteins from the lysate during the initial centrifugation of the affinity resin complexes.