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In Vivo Imaging of DiR Stained Spleen T-cell Distribution


Figure 1. DiR stock was prepared by dissolving 25 mg in 3 mL ethanol. Working solution of 320 $\mu\text{g}/\text{mL}$ was prepared by diluting 199 μL of stock solution in 5 mL PBS. T-cells isolated from the spleen were incubated with 320 $\mu\text{g}/\text{mL}$ DiR. After 30 min incubation, cells were spun down for 3 min at 1000 rpm at 4 $^{\circ}\text{C}$ resulting in a blue pellet. Cells were washed twice in PBS and injected intravenously (5×10^6 cells/mouse). Control group was injected with 5×10^6 cells/mouse in PBS. Mice were imaged with IVIS Spectrum at 10 min, 1hr, 6hr and 24 hrs post injection. Ideal filter set for DiR imaging is 710 nm excitation and 760 nm emission. Mice were imaged dorsally as well as ventrally at all time points. Brain, bones, spleen, liver, lungs and kidneys were harvested for ex vivo imaging 24 hrs post injection. Non-invasive in vivo imaging showed the homing process of injected T cells to the liver and spleen in real time, which was confirmed by ex vivo imaging.

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

Figure 2. T-cells isolated from the spleen were fluorescently stained with DiR and i.v. injected (5×10^6 cells/mouse) into a Nu/Nu mouse. Images above taken 24hrs post injection with IVIS Spectrum show cells homing to the spleen.

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Kalchenko et al., Use of lipophilic near-infrared dye in whole-body optical imaging of hematopoietic cell homing. *Journal of Biomedical optics*, September/October 2006, Vol 11(5).

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