

M5 Anti-HA Tag Rabbit Polyclonal Antibody

抗 HA 标签兔多克隆抗体使用说明书

产品名称	单位	货号
M5 Anti-HA Tag Rabbit Polyclonal Antibody	20 μ l	MF084-R-T
M5 Anti-HA Tag Rabbit Polyclonal Antibody	100 μ l	MF084-R-01

【STORAGE】

Store at -20°C. Stable for one year from the date of shipment.

【BACKGROUND】

该抗体高度特异识别重组蛋白 C 末端或 N 末端的 HA 标签，而不受邻近氨基酸组成的影响，可用于检测各种表达载体表达的 HA-Tag 融合蛋白。

【SOURCE】

This polyclonal antibody is produced by immunizing animals with a synthetic peptide containing the influenza hemagglutinin epitope (YPYDVPDYA) (KLH-coupled).

【SPECIFICITY】

HA-Tag Rabbit Polyclonal antibody detects over-expressed proteins containing the HA epitope tag.

【REACTIVITY】

All

【ISOTYPE】

Rabbit IgG1

【Concentration】

1mg/ml

RECOMMENDED ANTIBODY DILUTIONS

Western blotting	1:1000-5000
ELISA	1:1000-20000

For Western blots, incubate membrane with diluted antibody in 5% w/v nonfat dry milk, 1X TBS, 0.05% Tween-20 at 4°C with gentle shaking, overnight.

APPLICATION DATA

【备注】

本产品仅供科研使用。在确认产品质量出现问题时，本公司承诺为客户免费更换等量的质量合格产品。