

M5 Hiper Goat Anti-rat IgG/FITC antibody

使用说明书

产品名称	单位	货号
M5 Hiper Goat Anti-rat IgGFITC antibody	100ul	MF761-01

Storage: Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

产品介绍 Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10¹¹ variants.

抗体来源	Goat
克隆类型	Polyclonal
交叉反应	rat
产品应用	Flow-Cyt=1:100-1000 IF=1:100-1000, not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
分子量	150kDa
性状	Lyophilized or Liquid
浓度	1mg/ml
免疫原	Full length plasma protein
亚型	IgG
纯化方法	affinity purified by Protein A
储存液	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

【备注】

本产品仅供科研使用。在确认产品质量出现问题时，本公司承诺为客户免费更换等量的质量合格产品。