

M5 Goat anti-Chicken IgG whole serum antibody

Product	Unit	Cat.#
M5 Goat anti-Chicken IgG whole serum antibody	1ml	MF623-01

【Storage】 : Shipped at 4°C. Store at -20°C.

【产品介绍】 :

Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 8-17 mg/mL in adult blood. IgG is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.

- Goat Polyclonal whole antiserum , Unpurified.
- Reacts with: chk.
- Immunogen: Chicken IgG purified purified from serum.
- Predicted Molecular Weight: 180kDa.

【免疫原】 : Chicken IgG purified purified from serum

【缓冲液】 : 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol

【应用范围】 :

WB=1:10000-50000

ELISA=1:10000-100000

IHC-P=1:10000-50000

IHC-F=1:10000-50000

not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"